

*Anywhere you
want to go!*

hugo[®]

*Freedom, mobility
and confidence to
stay active and
go anywhere,
anytime!*

*the
rolling walker
with a seat*




*take your break
in style!*



*enjoy a free
and easy lifestyle!*





- A. Handlebar
- B. Backrest
- C. Brake Cables
- D. Seat
- E. Tie Rod
- F. Front Wheels

- G. Rear Wheels
- H. Folding Bracket
- J. Guide Pin
- K. Hand Screw
- L. Spring Pin
- M. Rear Leg Hinge
- N. Handgrip

SAVE THESE INSTRUCTIONS

To ensure your safety in using the HUGO® Rolling Walker, these instructions MUST be followed:

To assemble your HUGO® Rolling Walker, refer to the diagram on page 2 of this brochure and follow these step-by-step instructions prior to operating. For illustrated step-by-step instructions, please refer to the enclosed Quick Start Guide.

1. Remove the HUGO® Rolling Walker from the carton and examine all parts for shipping damage. DO NOT use the rolling walker if there is apparent damage or if the following parts are missing: **4 Triangular Hand Screws, 2 Front Wheels, 1 Backrest, 2 Handlebars.** If there are missing or damaged parts call HUGO® customer care, toll free, 1-866-321-HUGO (4846).
2. Lay the rolling walker down so the seat is facing the floor. Cut the plastic ties from the Rear Wheels (G) and unfold the legs. Straighten the legs until the Rear Leg Hinges (M) lock into place.
3. Flip the rolling walker over so that the seat faces up. Insert Front Wheel (F) into frame, orienting with the Guide Pin (J). Be sure to fully seat the Wheel until both Spring Pins (L) lock into place. Secure in place using the triangular Hand Screw (K). **HAND TIGHTEN ONLY.** Repeat for other Front Wheel.
4. Unfold HUGO® by pushing down on the Tie Rod (E) to spread apart the frame until the Folding Brackets (H) are in full extension. Push the Seat (D) down into place.
5. Insert the Handlebars (A) so that the Handgrips (N) are aligned with the Rear Wheels (G). The Brake Cables (C) should be on the **outside** of the Handlebars and Frame.
6. **ADJUST THE HANDLEBARS FOR PROPER FIT:** Stand upright behind the rolling walker with shoulders relaxed, arms hanging at your side, and eyes looking straight ahead. Standing in this position, the Handgrips (N) of the rolling walker should be adjusted to a height even with the wrists. Secure the Handlebars at the closest hole to the desired height. Insert and tighten the triangular Hand Screws (K) until secure. **HAND TIGHTEN ONLY.**
7. Insert the Backrest (B) so that the padded portion faces the front of the unit. Push down the buttons on the side of the backrest until the Spring Pins lock into place.

READ THESE IMPORTANT SAFETY INSTRUCTIONS PRIOR TO OPERATING YOUR HUGO® ROLLING WALKER:

1. Consult your doctor, therapist or other qualified professional to determine if the HUGO® Rolling Walker will suit your particular needs.
2. **IMPORTANT-** Before using, make sure that all parts are secure and the rolling walker frame is fully opened and in locked position.
3. To Use the Rolling Walker: For proper positioning and use of your HUGO® Rolling Walker, imagine a line drawn across the back of the rear wheels. Part of one foot should always be slightly ahead of that line as you walk.
4. To Operate the Brakes:
 - To slow down, squeeze the brake handles up.
 - To lock the brakes, push down on the brake handles until you feel or hear the brakes click into locked position. If the rear wheels of the rolling walker move while the brakes are in locked position, DO NOT use the rolling walker. Please call our toll-free number listed on page 4 for return or repair.
 - To release the brakes from locked position, simply pull up on the brake handles.
5. To Use as a Seat: Push your HUGO® Rolling Walker into position where you want to sit, then lock the brakes before sitting.

IMPORTANT- Before using the seat portion of your rolling walker, ALWAYS lock the brakes FIRST. DO NOT navigate the rolling walker while sitting on the seat. The brakes should ALWAYS be locked when the seat is being used. Risk of fall and serious injury may result if brakes are not locked while the user is sitting on the seat.

- DO NOT use the seat without the backrest. Falls are possible if the user leans or slides backward on the seat.



Safety Instructions *(continued)*

- NEVER attempt to sit on the seat when the rolling walker is on an incline. Risk of fall and serious injury may result.
 - Your HUGO® Rolling Walker is for individual use only and is NOT to be used as a wheelchair. DO NOT attempt to push the rolling walker while someone is sitting in the seat. Damage to the rolling walker frame or wheels, and serious injury to the user, could result from improper use.
 - **CAUTION-** DO NOT exceed the following weight limits, or serious injury could result:
Handles: 200 lbs.; Seat: 250 lbs.
 - The seat is not intended to transport people or objects. DO NOT use the seat to carry or move anything.
6. To Collapse for Storage or Transport: Pull up on the tie-bar until the frame is fully collapsed.

Your HUGO® Rolling Walker should be checked periodically to ensure the brakes are functioning properly and that all nuts and bolts are secure.



Important Information

Be sure to fill out the warranty registration card and mail it in. You must include your serial number, which can be found on the bar beneath the seat.

Keep your original purchase receipt in a safe place.

Limited Warranty

Your HUGO® Rolling Walker is warranted to be free from defects in materials and workmanship under normal use and service for 3 years on the frame and 1 year on all other parts. See enclosed Limited Warranty for complete warranty terms and conditions.

Contact us at 1-866-321-HUGO (1-866-321-4846), Monday-Friday, 8am - 6pm EST, if you encounter any problems or for information regarding repairs or maintenance.



Walking Tips

1. Before walking, allow plenty of time to perform a few simple stretching movements. As a result, you will be able to walk further and avoid injury.
2. Choose proper footwear. Wear walking shoes that are stable from side to side and well cushioned, snug in the heel, yet allow enough room to wiggle your toes freely.
3. Wear comfortable clothing. When walking near a road, it's wise to wear bright clothing that will be easily seen by passing motorists.
4. Bring a supply of water with you. It is important to stay hydrated not only on hot days, but at any time in dry climates.
5. Avoid walking in very hot or very cold weather. On those days, choose to walk at your local mall or indoor track.
6. Walk on soft ground. Softer ground is more foot-friendly than cement or asphalt. If possible, walk on grass or dirt paths that are flat, even and well manicured.
7. Exercise smart; know your limits. Consult your physician before beginning a fitness program, especially if you are substantially overweight, physically inactive and become easily fatigued, and / or have a history of heart disease, poor circulation or diabetes.

This product contains dry natural rubber.

Made in China

Distributed by
Access Product Marketing, LLC
Alpharetta, GA 30022

HUGO® is a trademark of
Access Product Marketing, LLC